



PRODUCT SPECIFICATION SHEET

WATER

HPLC GRADE

0.22 micron filtered

Main Catalog #: 23200HPLC-Package Size*

Available in the following sizes:

*Refer to Master Price List – Individual package sizes have unique size codes

Test	Specification	Typical Result
Fluorescence at emission maximum for Impurities at 450 nm	0.1 ppb max.	<0.1 ppb
LC Gradient Elution Suitability at 254 nm	To Pass Test	Pass
Residue after Evaporation	1 ppm max.	<0.1 ppm
Specific Conductance at 25C ¹	1.5 µS/cm max.	0.05 µS/cm
Optical Absorbance at: 190 nm	0.01 max.	0.00
200 nm	0.01 max.	0.00
210 nm	0.01 max.	0.00
250 nm	0.005 max.	0.001
400 nm	0.005 max.	0.001
Substances reducing Permanganate	To Pass Test	Pass

Form: Water. HPLC. Rev 1.1, 02/12, AM

¹ Meets stated value at the time of manufacture.

Not intended for use as a sterile or potable product.